



Support Surface Cover and Core: Working Together In Sweet Harmony

A study into how different types of polyurethane-coated covers impact the pressure redistribution and microclimate properties of medical support surfaces.



Aim

This study shows that not all fabrics are the same by demonstrating how properties of the support surface can be influenced when the fabric and core operate together.

BACKGROUND

Pressure redistribution has often been reported in relation to medical support surfaces. As far back as 1993, mattress covers were shown to have an impact¹. To date, however, the attention surrounding mattress performance has focused almost exclusively on the mattress contents, such as the properties of foam, the type of air system utilized, etc. with minimal attention being paid to the cover, except in regards to cleaning and care². Likewise, it is well documented that managing skin temperature and moisture levels is important for

maintaining skin integrity; for example, Fisher, Szymke and Apte (1978) report that a 1°C increase in skin temperature leads to a 13% increase in metabolic demand³. Relative humidity also affects the strength of the stratum corneum: at a relative humidity of 100% the stratum corneum is 25 times weaker than at 50% relative humidity⁴. However, there appears to be little research in bringing these factors together to establish the cumulative impact they have when considering the properties of a medical support surface as a whole.

Clinical Relevance

PRESSURE REDISTRIBUTION

Demonstrating how the mattress cover can have an impact on the redistribution of pressure, depending on the properties of the fabric¹.




IMMERSION & ENVELOPMENT

Demonstrating how the mattress cover fabric can have an impact on the extent of patient immersion and envelopment.









MICROCLIMATE CHARACTERISTICS

Demonstrating how different fabrics can influence the level of humidity in the area between skin and fabric.



Method

Three different tests were carried out on a range of support surface covers and cores in a laboratory setting: Pressure Mapping, Immersion & Envelopment and Temperature & Humidity


DARTEX® RANGE	PROPERTY OF COVER	DARTEX® RANGE	PROPERTY OF COVER
PRESSURE REDISTRIBUTION		MICROCLIMATE CHARACTERISTICS	
 CAR676	LOW STRETCH	 MIC200	HIGH BREATHABILITY 1000g/24hr/m ² , ASTM D1653 (g/m ² /24h)
 PER200	2 WAY STRETCH (High Stretch, Low Modulus)	 PER200	STANDARD BREATHABILITY Typical Polyurethane-coated fabric 500g/24hr/m ² , ASTM D1653 (g/m ² /24h)
 PER 406	4 WAY STRETCH (High Modulus)		
 END409	4 WAY STRETCH (High Stretch, High Modulus)		

FORCE/EXTENSION


The fabrics chosen have very different stretch properties, as can be seen on the Force/Extension graphs below:

These graphs demonstrate how much each fabric stretches when force is applied. On a support surface, the fabric needs to stretch to envelop the patient, and allow them to be immersed into the mattress.

FORCE/EXTENSION CURVES – WIDTH



FORCE/EXTENSION CURVES – LENGTH



For each mattress, two covers were made from each fabric: a tight cover and a loose cover.

Three different mattresses were specified for the experiment:

- Mattress 1** High specification foam mattress
- Mattress 2** Castellated high specification foam mattress
- Mattress 3** Foam and gel hybrid mattress

Covers for the mattresses were made using four fabrics from the Dartex Healthcare range:

- CAR676** Low stretch fabric
- PER200** 2 way stretch fabric (high stretch, low modulus)
- PER406** 4 way stretch (high modulus)
- END409** 4 way stretch (high stretch, high modulus)

Pressure Redistribution

PRESSURE MAPPING

A loose cover and a tight cover was made in each of the four fabric types. The covers were then applied to each of the three different mattresses and pressure mapped using an XSENSOR ForeSite SS Pressure Mapping system, with a 73kg adult male subject. The standard cover supplied with each mattress was also tested.



IMMERSION TEST

Using a system of weights and pulleys, a wooden mannequin weighing 178lbs was lowered onto the mattress. Percent immersion was calculated based on the mannequin thickness. The mannequin was raised from the support surface for a recovery of 300 seconds \pm 15 seconds between each trial. Three measurements were taken each time and the average recorded.



ENVELOPMENT TEST


The test uses a bulbous indenter, designed to replicate the shape of the buttocks. The indenter is embedded with pressure sensors, and the results are reported across four zones – **ZONE 1** shown here in dark blue, **ZONE 2** in light blue, **ZONE 3** in maroon and **ZONE 4** in orange.

surface for 60 seconds and pressure readings taken. The more even the pressure distribution across the four zones, the better the envelopment.

The test environment had an ambient temperature of 23 \pm 2 °C and relative humidity 50% \pm 5% which is specified in ISO 554–1976(E).

A load of 22 \pm 1lbs (representing the pelvic load of a 50th percentile male) was applied to the test support

Left lateral 40 Trochanter



Microclimate Characteristics

TEMPERATURE & HUMIDITY TEST

Using temperature and humidity sensors, the microclimate of a 73kg seated volunteer was monitored to see how the use of seating fabrics with different moisture vapour permeability (MVP) levels affected the temperature and humidity experienced. Sensors were placed underneath the thighs and buttocks of the volunteer, above and below the cover fabric. The cover was designed with MIC200 on one side (blue) and PER200 on the other (black). Data was collected using a Body View system supplied by Inside Climate. The test was run for 1 hour.




Results

PRESSURE MAPPING

Real differences were seen between the tested fabrics. The difference between best and worst fabrics was greater for the more “technical” mattress. The “best” combination was different for each mattress. For the low

stretch fabrics, loose covers gave lower peak pressure and lower average pressure. For the high-stretch fabrics, the results were very similar. The average pressure was largely independent of fabric.

WORST COVER - STANDARD




5 15 25 35 45

The same mattress core provided the best and worst results.

Coefficient of variation	162.22
Average (mmHg)	21.11
Peak (mmHg)	45.42
Area (cm ²)	3377.01

BEST COVER – PER 200 TIGHT



5 15 25 35 45


Coefficient of variation	130.47
Average (mmHg)	19.22
Peak (mmHg)	38.64
Area (cm ²)	3044.35

IMMERSION

The absolute level of **immersion** achieved is primarily determined by the mattress core, but for a given core, the maximum immersion can be increased by choosing the right cover fabric and the right tension i.e. tight or loose cover.

For each mattress, changing the cover fabric influences both the maximum pressures recorded and the distribution of pressure across the indenter. The stretch properties of the fabric have a large influence on the way the support surface is able to **envelop**:


KEY: CAR676 Tight/CAR676 Loose/PER200 Tight/PER200 Loose/PER409 Tight/PER409 Loose/END409 TIGHT/PER409 Loose



Sample	Mattress 1 % Immersion	Mattress 2 % Immersion	Mattress 3 % Immersion
CAR676 Tight	26.40	29.70	40.70
CAR676 Loose	26.80	29.20	40.70
PER200 Tight	28.70	30.50	41.90
PER200 Loose	29.30	30.20	48.00
PER406 Tight	26.00	31.10	41.50
PER406 Loose	28.40	29.90	41.00
END409 Tight	28.30	30.50	40.80
END409 Loose	27.00	28.90	38.00


ENVELOPMENT

MATTRESS 1 ZONE 1/ZONE 2/ZONE 3/ZONE 4




Cover	Zone 1 (mmHg)	Zone 2 (mmHg)	Zone 3 (mmHg)	Zone 4 (mmHg)
CAR676 Tight	80.39	62.45	10.08	0.16
CAR676 Loose	69.43	49.56	15.36	-0.21
PER200 Tight	53.57	51.80	7.19	-2.50
PER200 Loose	49.24	44.44	17.98	1.16
PER406 Tight	57.00	52.36	21.68	0.12
PER406 Loose	56.64	49.85	19.20	1.72
END409 Tight	63.58	55.61	17.23	1.70
END409 Loose	58.71	51.77	20.22	0.11

MATTRESS 2 ZONE 1/ZONE 2/ZONE 3/ZONE 4



Cover	Zone 1 (mmHg)	Zone 2 (mmHg)	Zone 3 (mmHg)	Zone 4 (mmHg)
CAR676 Tight	61.62	57.74	15.61	-0.02
CAR676 Loose	62.90	56.07	17.48	0.35
PER200 Tight	45.75	49.65	19.26	0.06
PER200 Loose	36.76	49.60	23.36	0.77
PER406 Tight	64.70	49.10	22.59	6.84
PER406 Loose	56.59	49.95	23.07	0.98
END406 Tight	65.19	57.08	21.14	0.10
END406 Loose	62.39	51.40	20.61	0.36

MATTRESS 3 ZONE 1/ZONE 2/ZONE 3/ZONE 4




Cover	Zone 1 (mmHg)	Zone 2 (mmHg)	Zone 3 (mmHg)	Zone 4 (mmHg)
CAR676 Tight	63.73	51.62	17.70	6.19
CAR676 Loose	48.92	48.89	14.37	6.21
PER200 Tight	45.58	44.47	14.75	7.73
PER200 Loose	55.60	46.66	17.69	11.26
PER406 Tight	59.31	49.02	16.95	10.02
PER406 Loose	33.91	48.31	22.20	9.57
END406 Tight	56.10	49.79	19.67	5.48
END406 Loose	39.75	47.20	20.00	7.07

MICROCLIMATE CHARACTERISTICS


The temperature profile was consistent, but the humidity levels in between skin and fabric (below) were significantly lower for the Microclimate fabric. When measuring absolute

humidity below the fabric, a greater increase in humidity was recorded for the Microclimate fabric; demonstrating the increased movement of moisture vapour.

ABSOLUTE HUMIDITY MIC200/PER200
(Above - between skin and cushion)



ABSOLUTE HUMIDITY MIC200/PER200
(Below - Under the cover)



Conclusion

The results infer that in the best medical devices, the support cover and core work harmoniously for the best patient outcome.

Moreover, the performance of the support surface core will depend on its compatibility with the cover. Where these separate parts are designed together, the support surface will perform to its full potential. Therefore, it is important that a like-for-like cover replacement (material, fit and moisture handling properties) takes place if needed during

its operational life to ensure the integrity of the support surface system is retained.. Choosing a fabric with high MVP can slow the humidity rise over time, keeping the patient feeling cooler and drier for longer. Further research is required to demonstrate the positive effects of the core and cover working together.



Trelleborg is a world leader in engineered polymer solutions that seal, damp and protect critical applications in demanding environments. Its innovative engineered solutions accelerate performance for customers in a sustainable way. The Trelleborg Group has local presence in over 40 countries around the world.



twitter.com/dartexcoatings
youtube.com/Trelleborg
linkedin.com/company/trelleborgecf
www.dartexcoatings.com

The research for this paper was undertaken by:

Richard Haxby, Katie Pearce, Ian Scott, Clare Tittershill, Tessa Turton & Claire Williams.

It has been presented at: European Seating Symposium in Dublin, Ireland (June 2016)

National Pressure Ulcer Advisory Panel (NPUAP) Annual Conference in St. Louis, MO, USA (March 2019)

References:

- 1 PS1 Foam Mattresses: A comparative evaluation special issue (1993) Medical Services Directorate.
- 2 British Healthcare Trades Association. (2012). BHTA guidance on the care, cleaning and inspection of healthcare mattresses. Available: <http://www.bhta.net/sites/default/files/document-upload/Protect%2C%20Rinse%2C%20Dry%20v2.pdf>. Last accessed 12th Aug 2016.
- 3 Fisher SV, Szymke TE, Apte SY, Kosiak M. Wheelchair cushion effect on skin temperature. Arch Phys Med Rehabil 1978; 59(2): 68-72.
- 4 Brienza DM, Geyer MJ. Using support surfaces to manage tissue integrity. Adv Skin Wound Care 2005; 18: 151-57.

The authors would like to acknowledge the following companies for the provision of equipment to undertake this study:

Direct Healthcare Services, EC Services Inside Climate GmbH, Stryker and Sumed.

© 2019 Trelleborg Group. At time of publication the information contained in this literature is believed to be correct and accurate.